

Hemostasis Innovation is Here. ►



Centralized management
for Hemostasis testing

Support laboratory decision-making

from any PC or tablet via secure network connection



“HemoHub Intelligent Data Manager is simple, intuitive, and modern—and was quickly adopted by the lab staff.”

— Dr. Inès Harzallah, Pharmacist Biologist, Head of Hemostasis Laboratory,
Group of Mulhouse and the South Alsace Region (GHRMSA), Mulhouse, France

Make the connection with **HemoHub Intelligent Data Manager.**

Enabling centralized information management, HemoHub connects across all Werfen Hemostasis testing systems, enhancing testing efficiency, streamlining workflows, and, ultimately, improving patient care.



Connectivity—standardize and optimize Hemostasis testing across systems and networks, for flexibility and efficiency, while supporting network strength and security

Quality control (QC) management—elevate quality management with comprehensive features, including Compare Analyzer QC and Compare Internal QC, and with AccuTrak 2.0™* peer-comparison review

Lab-decision support—access key insights, metrics, and key performance indicators (KPIs); flexible and robust rules support customer-defined and flexible-rules capabilities, performance tracking, and standardized processes

Management of HemoCell® Specialized Lab Automation—optimize workflows and improve diagnostic outcomes through centralized management

*AccuTrak2.0 is a peer-comparison portal. Users can easily download the QC files from HemoHub and upload them to AccuTrak 2.0.

Connectivity

Standardize Hemostasis testing across an entire network,
for greater flexibility and efficiency



Displays status of all connected instruments from a single dashboard

Enables monitoring of results, in real time, for entire network

Offers a user-friendly interface with color codes, alarms, and flags

Records all actions taken on critical results

Supports security and enables encrypted communication

QC management

Manage QC simply and efficiently,
from one central location



Displays and compares QC performance across multiple systems

Compares internal QC and reports generation (e.g., Bias, Uncertainty), enabling valuation of key metrics

Provides input data for AccuTrak 2.0

Allows performance verification studies†

†Only applicable for ACL TOP® Family 70 Series systems.

Lab-decision support

Helps enable rapid lab decisions through easy access to patient result information



Reduces manual review time

Displays, saves, and compares reaction curves via Customizable Curve Library

Provides preanalytical data to support decision-making


Allows rules configuration to enhance, standardize, and simplify operator workflows

Offers easy access to key performance indicators for performance monitoring

Management of HemoCell Specialized Lab Automation

Designed to address the unique needs of your laboratory

HemoCELL
SPECIALIZED
LAB AUTOMATION



HemoCell combines HemoHub, ACL TOP 770/750 LAS testing systems, and HemosIL reagents.

Centralized management through seamless integration


End-to-end automation, including automated performance verification tools

Intelligent sample routing and tracking with real-time adaptations

Customized sorting and sample banking

Enhanced performance, workflow, and lab-decision support

Standardizes and optimizes Hemostasis testing across systems and networks—for flexibility, efficiency, and improved patient care.



HIS/LIS/EMR = hospital information system; laboratory information system; electronic medical record

Corporate Headquarters

Plaza de Europa, 21-23
08908 L'Hospitalet de Llobregat
Barcelona, Spain
+34-93-4010101

Hemostasis and Acute Care Diagnostics

Headquarters and Technology Center

180 Hartwell Road
Bedford, MA 01730-2443
USA
+1-781-861-0710

Worldwide Locations

The Americas

Brazil
São Paulo
+55-11-46227878

Canada
Richmond Hill, ON
+1-800-552-2025 x6115

Colombia
Bogotá
+57-15-221-052

Mexico
Mexico City
+52-55-5262-1760

Uruguay
Montevideo
+5982-481-81-33

USA
Bedford, MA
+1-781-861-0710

Asia-Pacific

Australia
Sydney
+61-02-9098-0200

China
Beijing
+86-10-59756055

Hong Kong
+852-2792-7773

Shanghai
+86-21-66308671

India
New Delhi
+91-490-29-550

Japan
Tokyo
+81-3-5419-1301

Korea
Seoul
+82-1899-9217

Thailand
Bangkok
+66-271-226-28/9

Europe

Austria
Vienna
+43-1-256-58-000

Belgium
Brussels
+32-2-7252052

Czech Republic
Prague
+420-246-090-931

France
Paris
+33-182-30-86-00

Germany
Munich
+49-89-909070

Hungary
Budapest
+36-1-882-73-10

Italy
Milan
+39-02-25221

Lithuania

Kaunas
+370-37-313157

The Netherlands

Breda
+31-76-5480100

Poland

Warsaw
+48-22-336-18-00

Portugal

Lisbon
+351-214247312

Russia

Moscow
+7-499-124-45-59

Spain

Barcelona
+34-902-20-30-90

UK

Warrington, England
+44-1925-810141

werfen.com

For more information, contact your local Werfen sales representative.

ACL, ACL AcuStar, ACL Elite, ACL TOP, HemoCell, HemoHub, HemosIL, ReadiPlasTin, RecombiPlasTin, SynthAFax, and SynthASil are trademarks of Instrumentation Laboratory Company (d.b.a. Werfen) and/or one of its subsidiaries or parent companies and may be registered in the United States Patent and Trademark Office and in other jurisdictions. The Werfen logo is a trademark of Werfen and may be registered in the Patent and Trademark Offices of jurisdictions throughout the world. All other product names, company names, marks, logos, and symbols are trademarks of their respective owners.

©2025 Instrumentation Laboratory. All rights reserved.